L3-SERIES PIPETTE-L

MULTI MATE PIPETTES 8 & 12 CHANNEL

DEUTSCH	6 – 23
---------	--------

ENGLISH	24 – 39
---------	---------



FRANÇAIS	40 – 57
----------	---------



slamed[®]
instruments for sciences

L3-SERIES PIPETTE-L


MULTI MATE PIPETTES 8 & 12 CHANNEL


3

D



E

F





L3-SERIES PIPETTE-L

MULTI MATE PIPETTES 8 & 12 CHANNEL


5

5A
2 μ l - 1000 μ l

5C
5000 μ l, 10000 μ l

INHALT

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CONTENTS

1 - INTRODUCTION
2 - SETTING THE VOLUME
3 - ASPIRATING AND DISPENSING LIQUIDS
4 - PRE-RINSING
5 - DENSE AND VISCOUS LIQUIDS
6 - PIPETTE TIPS
7 - RECOMMENDATIONS
8 - RECALIBRATION
9 - TROUBLESHOOTING
10 - CLEANING AND STERILIZATION
11 - PIPETTE KIT
12 - SPARE PARTS

1 - INTRODUCTION

The pipette is a volumetric instrument designed to measure and transfer liquids precisely and safely. It can measure and transfer volumes from 0.1 μl to 10000 μl depending on the model.

The adjustable multichannel pipettes have been designed for the filling of laboratory microplates. The pipettes enable precise and simultaneous delivery of 8 or 12 preset-volume doses of liquid. The multichannel pipettes are produced in four ranges of volumes: 0.5-10 μl , 5-50 μl , 20-200 μl , 50-300 μl .

The pipettes are equipped with a digital counter which shows the pipetting volume. The volume setting is done by turning the pipetting pushbutton knob (Fig. 1A2) or the black adjustment knob (Fig. 1B) in the right direction. The volume range is shown on the pipetting pushbutton (Fig. 1A1).

The volume range from 0.1 μl to 10000 μl is covered by 10 pipette models.

Pipette volume	Application
L2, L10	Measurement and transfer of micro-volumes, DNA sequencing and enzyme-assay applications.
L20, L50, L100, L200, L250, L300, L1000	Measurement and transfer of general aqueous solution, acids and bases
L5000, L10000	Measurement and transfer of large volumes.

The liquid is aspirated into disposable tips attached to the pipette shaft (Fig. 1C). Disposable tips ensure maximum safety and eliminate possibility of cross-contamination of the liquid samples.

To protect the user from contaminated tips the pipette is equipped with a built-in tip ejector, (Fig.1D).

Single channel pipette's ejector

The construction of the ejector enables the user to set up the length. The adjustable tip ejector accommodates every variety of tips available on the market. When using narrow tubes, it may be necessary to remove the tip ejector. It is simply removed by pulling down.

Adjusting of the Tip Ejector Length (2-1000 μ l)

(Fig. 5A).

Adjust the tip ejector length by turning the ejector cap after removing the ejector pushbutton. To lengthen the tip ejector turn the ejector cap clockwise. To shorten the tip ejector length turn the ejector cap counterclockwise.

Adjusting of the Tip Ejector Length 5000 and 10000 μ l

(Fig. 5B).

The length of tip ejector is adjusted by screwing in or screwing out the tip ejector stem with a screwdriver. Turn the screwdriver counter-clockwise to increase the length of tip ejector, and clockwise to reduce the length of tip ejector. The ejector can be adjusted in the range of 5 mm.

If above described method of ejector adjustment is not sufficient or the diameter of the ejector opening is too large to eject the tip it is necessary to put the ejector cap "M" onto the ejector, (Fig. 5C).

Multichannel pipette's ejector

The multichannel pipette construction allows to rotate the multichannel module and to eject the tips through 360°, therefore, the most convenient position may be selected to fill the microplates. It is recommended to turn the module

L3-SERIES PIPETTE-L

MULTI MATE PIPETTES 8 & 12 CHANNEL

clockwise. Due to long stroke of ejector in multichannel pipette, the adjustment of its length is not necessary.

The pipette is a high quality instrument which offers excellent accuracy and precision.

The accuracy and precision (repeatability) of liquid volume depend on the quality of tips used. The values for accuracy and precision given in the table below were obtained using high quality tips.

The single channel pipettes:

Model	Cat. no.	Volume [μ l]	Accuracy [%]	Precision [%]	Tip μ l
L2	47 70 001	0.2	± 12.0	± 6.0	10 μ l
		1.0	± 2.7	± 1.3	
		Max 2.0	± 1.5	± 0.7	
L10	47 70 002	Min 0.5	± 4.0	± 2.8	10 μ l
		5.0	± 1.0	± 0.6	
		Max 10.0	± 0.5	± 0.4	
L20	47 70 003	Min 2	± 3.0	± 1.5	200 μ l
		10	± 1.0	± 0.5	
		Max 20	± 0.8	± 0.3	
L50	47 70 004	Min 5	± 2.5	± 2.0	200 μ l
		25	± 1.0	± 0.6	
		Max 50	± 0.8	± 0.4	
L100	47 70 005	Min 10	± 1.6	± 0.80	200 μ l
		50	± 0.8	± 0.24	
		Max 100	± 0.8	± 0.20	
L200	47 70 006	Min 20	± 1.2	± 0.60	200 μ l
		100	± 0.8	± 0.25	
		Max 200	± 0.6	± 0.20	
L250	47 70 007	Min 50	± 1.0	± 0.4	250 μ l
		125	± 0.8	± 0.3	
		Max 250	± 0.6	± 0.3	
L1000	47 70 008	Min 100	± 1.6	± 0.40	1000 μ l
		500	± 0.7	± 0.20	
		Max 1000	± 0.6	± 0.15	
L5000	47 70 009	Min 500	± 1.2	± 0.50	5000 μ l
		2500	± 0.6	± 0.20	
		Max 5000	± 0.5	± 0.15	
L10000	47 70 010	Min 1000	± 2.5	± 0.6	10000 μ l
		5000	± 0.8	± 0.3	
		Max 10000	± 0.5	± 0.2	

The multichannel pipettes:

Model	Cat. no.	Volume [μ l]	Accuracy [%]	Precision [%]	Tip μ l
ME 8-10 MT 12-10	47 70 064 47 70 068	Min 0.5 5 Max 10	± 10.0 ± 4.0 ± 2.0	± 8.0 ± 2.0 ± 1.2	10
ME 8-50 MT 12-50	47 70 065 47 70 069	Min 5 25 Max 50	± 4.0 ± 3.0 ± 1.6	± 2.5 ± 1.2 ± 0.6	200
ME 8-200 MT 12-200	47 70 066 47 70 070	Min 20 100 Max 200	± 3.0 ± 1.5 ± 1.0	± 1.5 ± 0.8 ± 0.6	
ME 8-300 MT 12-300	47 70 067 47 70 071	Min 50 150 Max 300	± 1.6 ± 1.2 ± 1.0	± 1.5 ± 1.0 ± 0.6	300

The accuracy and precision are obtained, using a gravimetric method, performing at least 10 measurements of distilled water at a temperature of $20 \pm 1^\circ\text{C}$, according to EN ISO 8655 standard.

The pipette design enables the user to perform the recalibration process according to the rules presented in section 8.

2 - SETTING THE VOLUME

The counter displays three figures to be read from top to bottom. Additional to the figures on the lower wheel are printed graduations to enable a volume setting in the range increment of each pipette model.

Examples of the meaning of the black and red figures:

The single channel pipettes:

Model	Counter readings	Set volume	Increment
L2	1 2 5	1.25 μ l	0.002 μ l
L10	0 7 5	7.5 μ l	0.02 μ l
L20	1 2 5	12.5 μ l	0.02 μ l
L50	0 3 8	38 μ l	0.2 μ l

L3-SERIES PIPETTE-L

MULTI MATE PIPETTES 8 & 12 CHANNEL

Model	Counter readings	Set volume	Increment
L100	0 7 5	75 μ l	0.2 μ l
L200	1 2 5	125 μ l	0.2 μ l
L250	2 4 5	245 μ l	0.2 μ l
L1000	0 7 5	0.75 ml	2 μ l
L5000	1 2 5	1.25 ml	10 μ l
L10000	0 7 5	7.5 ml	20 μ l

The multichannel pipettes:

Model	Counter readings	Set volume	Increment
ME 8-10 MT 12-10	0 3 5	3.5 μ l	0.02 μ l
ME 8-50 MT 12-50	0 6 5	6.5 μ l	0.10 μ l
ME 8-200 MT 12-200	0 8 5	85.0 μ l	0.20 μ l
ME 8-300 MT 12-300	2 5 0	250.0 μ l	1.0 μ l

The volume of the pipette is set by the knob in the pipetting pushbutton (Fig. 1A2) or by the black adjustment knob (Fig. 1B). To attain the maximum accuracy, set volume must be approached from a higher value by diminishing counter readings.

If the desired value is lower than the previous value, simply adjust the value down to 1/3 turn above the required setting. If the desired value is higher than the previous

value, adjust the value to 1/3 turn above the required value. The setting should then be performed in the following manner:

Slowly turn the pipetting pushbutton knob or the black adjustment knob to the required volume.

If the knob is accidentally turned too far, below the required value, the process must be repeated.

3 - ASPIRATING AND DISPENSING LIQUIDS

Place a tip on the shaft of the pipette. See section 5 for the appropriate tip. Press the tip on firmly using a slight twisting motion to ensure a positive, airtight seal.

Important: Never aspirate liquids into the pipette without a tip attached.

2, 10	≤ 1
20, 50, 100	$2 \div 3$
200, 250, 300, 1000	$2 \div 3$
5000	$3 \div 6$
10000	$5 \div 7$

Aspiration

Press the pushbutton to the first positive stop, (Fig. 2A). Holding the pipette vertically, immerse the tip into the sample liquid. The depth to which the tip is immersed in the sample liquid depends on the model.

Release the pushbutton slowly and smoothly to aspirate the sample, (Fig. 2B). Wait one second and then withdraw the tip from the liquid. Wipe any droplets away from the outside of the tip using a medical wipe.

Avoid touching the orifice of the tip.

Dispensing

- Place the end of the tip against the inside wall of the vessel at an angle of 10 to 40 degrees.
- Press the pushbutton smoothly to the first stop, (Fig. 2C). Wait one second.
- Press the pushbutton to the second stop to expel any remaining liquid, (Fig. 2D).
- Keeping the pushbutton depressed to the very end, remove the pipette by drawing the tip against the inside surface of the receiving vessel.
- Release the pushbutton to its starting position, (Fig. 2E).

L3-SERIES PIPETTE-L

MULTI MATE PIPETTES 8 & 12 CHANNEL

- Eject the tip by pressing the tip ejector button, (Fig. 2F). Remember to change the tip whenever a different kind of liquid is to be sampled.

Filters

A replaceable filter installed in a seat in the bottom part of the shaft is used in 5000 μ l and 10000 μ l pipettes (Fig. 3L). The filter prevents the aspirated liquid from entering into the shaft and thus from polluting the inside of the shaft and the piston. Using the filter is especially important when aspirating and dispensing large volumes of liquid.

If the filter becomes wet during liquid aspiration it should be replaced with a new one.

4 - PRE-RINSING

When pipetting liquids of higher viscosity or lower surface tension than water (e.g. sera or organic solvents), a film of liquid is formed on the inside wall of the pipette tip. This film can create an error. Since the film remains relatively constant in successive pipetting operations with the same tip, this error can be avoided by forming the film before transferring the first sample. This is done by aspirating a sample and dispensing it back into the same vessel. Since the film is already formed, all of the following samples will have better accuracy and repeatability.

This pre-rinsing operation should be repeated when the volume to be aspirated is changed or when a new tip is used.

5 - DENSE AND VISCOUS LIQUIDS

The pipette specifications of accuracy and precision are based on pipetting distilled water. The handling of liquids with physical qualities of density, viscosity and surface tension differing extremely from water may need a gravimetrically checked compensation of the volume setting. Normally the degree of error resulting from heavy or viscous liquids is negligible if the pipetting is done slowly and carefully. It is most important to give the liquids some time to follow the change of pressure by holding the pipette tip in its position for at least 2 sec. after the aspiration and the blow out stroke.

If in extreme cases this method of operation does not result in accurate values, a compensation could be achieved as follows:

Weigh the liquid pipetted when the pipette is set to the nominal value. Then calculate the set-off from the nominal value:

$$\text{Corr: val.} = 2 \times \text{nom. val.} - \frac{m}{\gamma}$$

m - *weight of the sample*

γ - *density of liquid*

Check this operation once again and correct if necessary. Note the corrected value for further pipetting the same kind of liquid.

6 - PIPETTE TIPS

Ratiolab tips are made from high performance polypropylene and their quality guarantees the precision and accuracy associated with the pipette. Strict control is maintained throughout the manufacturing process to ensure the highest quality.

The use of inferior quality tips will seriously degrade the performance of the pipette.

Tips 10

These tips are used for volumes between 0.1 μl and 10 μl . They are used with the L2, L10, ME-10 and MT-10 models.

Tips 200

These tips are used for volumes between 2 μl and 200 μl . They are used with the L20, L50, L100, L200, ME-50, MT-50, ME-200 and MT-200 models.

Tips 250

These tips are used for volumes between 50 μl and 250 μl . They are used with the L250 model.

Tips 300

These tips are used for volumes between 50 μl and 300 μl . They are used with the ME-300 and MT-300 models.

Tips 1000

These tips are used for volumes between 100 μl and 1000 μl . They are used with the L1000 models.

Tips 5000

These tips are used for volumes between 1000 μl and 5000 μl . They are used with the L5000 models.

Tips 10000

These tips are used for volumes between 1000 μl and 10000 μl . They are used with the L10000 models, which are equipped with a white pushbutton.

7 - RECOMMENDATIONS

Observing the following recommendations will ensure maximal possible accuracy and precision of liquid sampling.

- Make sure to operate the pipette slowly and smoothly.
- The depth of immersion in the sample liquid should be the minimum necessary and should remain constant during aspiration.
- The pipette should be held in a vertical position.
- Change the tip when volume setting is changed or when a different liquid is to be aspirated.
- Change the tip if a droplet remains on the end of the tip from the previous pipetting operation.
- Each new tip should be pre-rinsed with the liquid to be pipetted.
- Liquid should never enter the pipette shaft. To prevent this:
 - Press and release the pushbutton slowly and smoothly.
 - Never turn the pipette upside down.
 - Never lay the pipette on its side when there is liquid in the tip.
- Never force the volume setting beyond its recommended limits.
- When pipetting liquids with temperatures different from the ambient temperature, it is recommended to pre-rinse the tip several times before use.
- Do not pipette liquids with temperatures above 70°C.
- When pipetting acids or corrosive solutions which emit vapours, it is recommended to disassemble the shaft and to rinse the piston and seal with distilled water after finishing the pipetting operation.

8 - RECALIBRATION

Pipettes are calibrated by gravimetric method, using high quality tips and distilled water, at the temperature $20 \pm 1^\circ\text{C}$, according to EN ISO 8655 standard.

If during pipette operation you find that the accuracy error (the difference between the real aspirated volume and the preset volume) exceeds the permissible value given in the table in section 1, the pipette recalibration procedure should be carried out.

Before starting the recalibration it is necessary to check

whether the following requirements have been fulfilled during error determination:

- the ambient temperature, and the temperature of the pipette, tips and water was identical
- the density of the liquid used is close to that of distilled water
- the balance with appropriate sensitivity has been used.

Volume checked [μl]	Balance sensitivity [mg]
0,1 - 10	$\leq 0,001$
10 - 100	$\leq 0,01$
> 100	$\leq 0,1$

- mg/ μl conversion factor has been taken into account
- the requirements given in sections 3 and 7 have been fulfilled

If the above conditions are satisfied and the accuracy error for selected volume given in section 1 exceeds the permissible value, the pipette recalibration procedure should be carried out.

The recalibration can be performed within one full turn of the key to the right or to the left only.

Recalibration conditions:

- Ambient temperature and the temperature of the pipette, tips and liquid should be within the range 20-25°C and stabilised during weighing within $\pm 0.5^\circ\text{C}$
- Measurements should be conducted using distilled water
- Balance sensitivity should be suitable for the volume to be controlled

Recalibration procedure:

- Set the dose volume depending on the pipette volume according to the following table:

L3-SERIES PIPETTE-L

MULTI MATE PIPETTES 8 & 12 CHANNEL

The single channel pipettes:

Model	Range of the pipette volumes [μ l]	Preset volume [μ l]	Permissible volumes [μ l]	Volume change ΔV for full turn of the calibration key [μ l] (24 increments)
L2	0.1 - 2	0.2	0.176 - 0.224	0.06
L10	0.5 - 10	0.5	0.48 - 0.52	0.33
L20	2 - 20	2	1.94 - 2.06	0.63
L50	5 - 50	5	4.875 - 5.125	2.50
L100	10 - 100	10	9.84 - 10.16	2.50
L200	20 - 200	20	19.76 - 20.24	6.30
L250	50 - 250	50	49.5 - 50.5	6.30
L1000	100 - 1000	100	98.4 - 101.6	25.00
L5000	1000 - 5000	1000	994 - 1006	125.00
L10000	1000 - 10000	1000	975 - 1025	250.00

The multichannel pipettes:

Model	Range of the pipette volumes [μ l]	Preset volume [μ l]	Permissible volumes [μ l]	Volume change ΔV for full turn of the calibration key [μ l] (24 increments)
ME 8-10 MT 12-10	0.5 - 10	0.5	0.45 - 0.55	0.33
ME 8-50 MT 12-50	5 - 50	5	4.8 - 5.2	1.67
ME 8-200 MT 12-200	20 - 200	20	19.6 - 20.4	6.30
ME 8-300 MT 12-300	50 - 300	50	49.2 - 50.8	10.00

- Perform 5 aspirations, weigh each one and calculate the average value of the aspirations.
- Perform three aspiration series (each series should include the aspirations from all channels), weigh each time and calculate the average value of the aspirations.
- Calculate average aspirated volume in [μ l] multiplying the average aspiration amount [mg] by the distilled water density coefficient [μ l/mg], which depends on temperature and pressure according to the following table:

Temperature [°C]	Pressure [kPa]		
	95.0	101.3	105.0
20	1.0028	1.0029	1.0029
21	1.0030	1.0031	1.0031
22	1.0032	1.0033	1.0033
23	1.0034	1.0035	1.0036
24	1.0037	1.0038	1.0038
25	1.0039	1.0040	1.0040

If the average aspirated volume exceeds the permissible value, the following should be done:

- Remove the pipetting pushbutton, (Fig. 4A),
Warning: The pipetting pushbutton consists of 2 parts: the knob (Fig. 1A2) and the pushbutton (Fig. 1A1). After removal of the pushbutton, both parts are separated.
- Holding the volume setting knob to protect it against rotation, insert the calibration key into the cuts of the calibration screw, (Fig. 4B),
- Turn the key clockwise to reduce the aspirated volume, or counter-clockwise to increase the volume, (Fig. 4C). One full turn of the calibration key changes the pipette aspiration volume by the amount given in the table.
- Take out the key and fix the pipetting pushbutton (Fig. 4D). The pushbutton should be fixed by placing first the knob on the arbor (Fig. 1A2) and then the pushbutton (Fig. 1A1).

Determine the average aspirated volume. The average volume should be within the permissible range given in the table. If the volume exceeds the values stated, the recalibration procedure should be repeated.

In case of pipetting the liquids with physical properties considerably different from those of water, follow the rules given in section 5.

9 - TROUBLESHOOTING

If the pipette aspirates too little liquid or there is a leak from the tip, check for the following possible causes and perform appropriate correcting measures:

L3-SERIES PIPETTE-L

MULTI MATE PIPETTES 8 & 12 CHANNEL

Problem	Cause	Solution
Droplets of liquid remain in the pipette tip.	The pipette tip is emptied too fast.	Decrease the speed of pressing the pipette pushbutton.
	The tip wettability has increased due to extensive use.	Replace the tip with a new one.
Bubbles of air appear in the liquid aspirated into the tip.	The pipette tip immersion is too shallow.	Immerse the tip deeper according to the instructions.
	The pipette tip is incorrectly pressed onto the pipette shaft.	Press the pipette tip firmly.
	The tip is damaged or worn out due to extensive use.	Replace the tip with a new one.
The pipette incorrectly aspirates the liquid or liquid drops out from the tip.	The pipette tip is incorrectly pressed onto the pipette shaft.	Press the pipette tip firmly.
	The shaft nut is loose (Fig. 3F) in the models 2 - 1000	Tighten the shaft nut.
	The sealing surface of the shaft is cracked or scored.	Remove the tip ejector. Unscrew the shaft nut, inspect the shaft and the piston assembly. Replace the damaged parts (see Section 12). When reassembling the pipette, the nut should be hand tightened. In the models L2, L10 and L20, the damage of the shaft may also cause a damage of the piston assembly. Replace the damaged parts (see Section 12). When reassembling the pipette, the nut should be hand tightened. To remove the tip ejector in models L5000 and L10000, remove the ejector pushbutton (Fig. 3N) and using a screwdriver unscrew the tip ejector by turning the screwdriver counter-clockwise.

Problem	Cause	Solution
The pipette incorrectly aspirates the liquid or liquid drops out from the tip.	Damage to the piston or seal due to prolonged use with the aggressive liquids.	Disassemble the pipette as described above. Replace the piston, seal and O-ring (see Section 12). Rinse the inside of the shaft in distilled water and dry. Lubricate the seal and O-ring with the lubricant, that is included with each pipette. Following piston replacement, a calibration procedure is required. Note: The parts of L2 and L10 pipette should be lubricated evenly with a minimum amount of lubricant.
	The pipette is re-assembled improperly.	Disassemble the pipette and reassemble it, observing the proper sequence of steps (Fig. 3).
	No lubricant on the sealing elements.	Remove the tip ejector Unscrew the shaft nut, remove the shaft, piston assembly, seal and O-ring. Rinse the removed parts in distilled water and dry thoroughly. Lightly lubricate the inside surfaces of the seal and the O-ring with the included lubricant. Reassemble the pipette in the reverse order.
	Contamination of the inside of the pipette caused by extensive aspiration of chemically aggressive liquids or because liquid got inside the pipette.	Remove the tip ejector. Unscrew the nut, remove the shaft, piston assembly, seal and O-ring. Rinse the removed parts with distilled water and dry thoroughly. Lightly lubricate the inside surfaces of the seal and the O-ring with the lubricant. Reassemble the pipette in the reverse order.

The shaft of the 5000 and 10000 models should be autoclaved without the filter.

The handle and the outside of the pipette may be cleaned using a cloth damped in isopropanol.

If the pipette malfunction persists after carrying out the above steps, send the pipette to your **Ratiolab** service representative.

Before returning the pipette, make sure that it is free from any chemical, radioactive or microbiological contamination

that might pose any danger during transportation and repairing. Clean the pipette as thoroughly as possible.

10 - CLEANING AND STERILIZATION

Cleaning

External surfaces of the pipetting pushbutton, the ejector pushbutton, the handgrip, the shaft nut and the adjustment knob may be cleaned using a cloth dampened in isopropyl alcohol. The remaining parts removed from the pipette during pipette disassembly may be washed with distilled water or isopropyl alcohol.

Sterilization:

The pipette can be sterilized in the autoclave at 121°C for 20 minutes. After sterilization, the pipette should be dried and cooled to room temperature.

It is recommended:

- to sterilize the pipettes in autoclave with an initial vacuum and drying cycle,
- prior to sterilization unscrew the shaft nut slightly in the L2-L1000 pipettes, and unscrew the shaft slightly in the L5000 and L10000. After autoclaving these parts should be screwed tight again.

Note: Pipettes should be autoclaved with the filter removed (L5000 and L10000). After autoclaving, check that the shaft nut is tightly fitted to the pipette handle

- prior to sterilization the multichannel pipette unscrew the nut connecting the handle and module slightly. After autoclaving these parts should be screwed tight again.

The precision of the results should not alter if the pipetting process and autoclaving are carried out as described in this manual. Because a slight change in the accuracy of the dosage may occur, it is recommended to:

- check the calibration of the pipette after the initial first, third and fifth autoclaving cycles and then after every 10 autoclaving cycles.

Ultra Violet (UV) Sterilization

The pipettes are UV resistant. However, the recommended distance from the radiation source to the exposed element should be not less than 50 cm. Prolonged and very intense UV exposure can cause decoloration of pipette parts without effecting their performance.

11 - PIPETTE KIT

The pipettes are delivered in the kits including:

- Pipette
- Instruction manual
- Calibration key
- Pipette holder
- Ejector regulation spacers
(for pipette models L2 – L1000)
- Ejector cap (for pipette models L5000, L10000)
- Identification labels
- Filters (for pipette models L5000, L10000)
- Lubricant

The pipette holder assembly diagram is shown in Fig. 5.

12 - SPARE PARTS

All the spare parts indicated in Fig. 1, 3, 4, 5:

A: Pipetting pushbutton A1: Pushbutton A2: Knob

B: Calibration knob

C: Shaft

D: Ejector

F: Shaft nut

G: Piston assembly

H: Spacer

I: O-ring

J: Seal

K: Calibration key

L: Filter

M: Ejector cap

can be ordered from your **Ratiolab** representative (type of pipette and name of the part for this pipette should be specified).

Warning: The replacement of the piston requires conducting of calibration procedure according to section 8.

All rights reserved.

Product described in this manual is subject technical modification. Measuring tolerances correspond to EN ISO 8655.

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