

PCR-Platten Kompatibilität / PCR plate compatibility

Ratiolab Bestell-Nr. / Order-no.		8610011	8610012	8610013	8610014	8610016	8610017	8610020	8610021
Ausführung / Description		96-Well 0.2ml fully skirted	96-Well 0.2 mL semiskirted frosted	96-Well 0.1 mL sub semi- skirted frosted	96-Well 0.2 mL nonskirted frosted	96-Well 0.2 mL sub semi- skirted frosted	96-Well 0.1 mL semiskirted white	384-well 40 µL fully skirted frosted	384-well 40 µL fully skirted white
Standard Well (0.2 mL height: 20.58 mm)			X		X	X			
Low Profile Well (0.1 mL height: 15.5 mm)				X			X		
Device	Note/ Holder No.								
Applied Biosystems									
7000			X			X			
7300	4367032		X			X			
7500	4367032		X			X			
7900HT Fast (with tray) (0.1mL type)				X					
7900HT Fast (with tray) (0.2mL type)			X			X			
StepOne	48-Well								
StepOne Plus				X			X		
7500Fast	4359652			X					
Via7 Fast (0.1mL model)	4459846			X					
Via7 Fast (0.2mL model)	4459845		X			X			
QuantStudio Fast 96-well Block (0.1mL model)	4459846			X					
QuantStudio 3/5/6/7 Fast 96-well Block (0.1mL model)	4379983			X			X		
QuantStudio Fast 96-well Block (0.2mL model)	4459845		X			X			
QuantStudio 3/5 Fast 96-well Block (0.2mL model)	4381850		X			X			
QuantStudio 6/7 (0.1mL model) Real-Time PCR System	384-Well							X	X
5700	Sequencer		X			X			
7700	Sequencer		X			X			
310/3100/3130	Sequencer		X			X			
3700/3730	Sequencer		X			X			
3730x/3500xI	384-Well							X	X
Sequencing Plates	384-Well								
3730, 3730xI, 3500, 3500xI	384-Well								
ProFlex™ PCR-System 384	384-Well								
Roche									
LightCycler 480							X		
LightCycler 480 II							X		
LightCycler 96							X		
LightCycler Nano	32-Well								
LightCycler 480 Multiwell Plate 384	384-Well								
Bio-Rad									
iCycler iQ		X			X				
iQ5		X			X				
MyiQ		X			X				
Chromo4		X			X				
Opticon I		X					X		
Opticon II		X					X		
mini Opticon	48-Well								
CFX96		X					X		
CFX384								X	X
Stratagene									
Mx4000/3000P/3005P					X				
Eppendorf									
MasterCycle ep Realplex					X				
Gradient Cycler ep		X						X	X
Thermo Scientific									
PikoReal 24 Real-Time PCR	24-Well								
Biometra									
Topcal Thermocycler					X				
qTOWER ³ B4 G 96 well	96 well				X				
qTOWER ³ B4 G 384	384-Well							X	X
Amersham									
Megabase 4000								X	
Qiagen®									
Rotogene®								X	X

Die obige Kompatibilitätstabelle dient nur zur Orientierung. Ratiolab haftet nicht für technische Änderungen oder Änderungen an den aufgeführten Geräten.
The compatibility chart above is given for guidance. Ratiolab is not responsible for any technical modifications or changes made to the devices listed.

Die Spezifikationen (= Präzision und Genauigkeit) sind nicht getestet.
The specifications (=inaccuracy and imprecision values) are not tested.