

 <small>disposables for sciences</small>	Qualitätsmanagement Handbuch	Formular P Form P	Seite page Seite 1 von 1
	Produkt-Spezifikation Product Specification	Änd. Datum Date of change 2022-12-05	Änd. Stand Revision 4

Allgemeine Produktdaten

ratiolab® Autoklavierbeutel - BIOHAZARD rot

Warnhinweis mit Indikationsfeld. Markante rote Beutel zur leichten Identifizierung und Entsorgung von kontaminiertem Material. Das weiße Indikationsfeld verfärbt sich durch den Autoklavierprozess dunkel. Mit Biohazard-Symbol und Sicherheitsaufdruck in Deutsch, Englisch, Französisch. Autoklavierbar bei 134°C, zum Autoklavieren darf der Beutel nur leicht verschlossen werden. Markierung maximale Befüllung (ausgenommen 700 32 30).

General Product Information

ratiolab® Autoclavable bags - BIOHAZARD red

Warning Instruction and sterilization indicator patch. Distinctive red autoclavable bags for easy identification, for the disposal of contaminated labware. Rectangle white indicator patch field turns into dark when autoclaved. With Biohazard symbol and safety imprint in German, English, French language. Autoclavable at 134°C, for autoclaving, the bag must not be completely closed. Marking maximum filling (except 700 32 30).

Material / Material: Polypropylen / Polypropylene

Warentarif-Nummer / Intrastat Number- HS Code: 39232990

Sterilisiert / pre-sterilized: Ja/Yes Nein/No

CE-Zeichnung / CE-marking: Ja/Yes Nein/No

B x H W x H mm	Toleranz tolerance	Volumen Volume ca.	Folienstärke Film Thickness	Packungseinheit Packaging	Bestell Nr. Order No.
300 x 500	+/- 5%	6 L	50 µm	1 x 500	70 03 200
400 x 780	+/- 5%	22 L	50 µm	1 x 500	70 03 210
600 x 800	+/- 5%	32 L	50 µm	1 x 500	70 03 220
600 x 800	+/- 5%	32 L	50 µm	1 x 100	70 03 225
700 x 1100	+/- 5%	72 L	50 µm	1 x 350	70 03 230
700 x 1100	+/- 5%	72 L	50 µm	1 x 75	70 03 235

This document was sent online and is valid without signature

Erstellt durch / Form edited by:	Aleksandra Mil / Auszubildende	Freigegeben von / Form approved by:	H. Bauer / QS/Technik
-------------------------------------	--------------------------------	--	-----------------------