

 <small>disposables for sciences</small>	Qualitätsmanagement Handbuch	Formular P Form P	Seite page Seite 1 von 1
	Produkt-Spezifikation Product Specification	Änd. Datum Date of change 2022-12-05	Änd. Stand Revision 4

Allgemeine Produktdaten

ratiolab® Kryo-Boxen, Format 133 x 133 mm

Aus tiefkältebeständigem Karton. Temperaturbeständig bis -80°C. Wahlweise wasserabweisend oder mit wasserfester Kunststoff-Spezialbeschichtung. Höchste Flexibilität durch freie Kombinationsmöglichkeiten mit verschiedenen Rastereinsätzen. Kryo-Boxen sind ohne Rastereinsatz auch für größere Gefäße verwendbar.

General Product Information

ratiolab® Cryo Boxes, Format 133 x 133 mm

Ultra low temperature-resistant cardboard. Temperature-resistant down to -80°C. Available with a water-repellent finish or with a waterproof plastic coating. Highest flexibility: Cryo Boxes can be combined with different Grid Inserts. The Cryo Boxes can be used for larger vessels when grid is removed.

Material / Material: Karton / cardboard

Warentarif-Nummer / Intrastat Number- HS Code 48192000

Sterilisiert / pre-sterilized: Ja/Yes Nein/No

CE-Zeichnung / CE-marking: Ja/Yes Nein/No

Farbe Color	Verpackungs- einheit Quantity per Pack	Höhe / Height 32 mm	Höhe / Height 50 mm	Höhe / Height 75 mm	Höhe / Height 100 mm	Höhe / Height 130 mm
		Bestell-Nr. Order No.	Bestell-Nr. Order No.	Bestell-Nr. Order No.	Bestell-Nr. Order No.	Bestell-Nr. Order No.
mit wasserabweisender Standardbeschichtung / Water-repellant standard coating						
weiß / white	10	51 31 321	51 31 501	51 31 751	51 31 001	51 31 131
blau / blue	10	51 31 322	51 31 502	51 31 752	51 31 002	51 31 132
rot / red	10	51 31 323	51 31 503	51 31 753	51 31 003	51 31 133
grün / green	10	51 31 324	51 31 504	51 31 754	51 31 004	51 31 134
gelb / yellow	10	51 31 325	51 31 505	51 31 755	51 31 005	51 31 135
mit Wasserfester Kunststoffbeschichtung / with waterproof plastic coating						
weiß / white	10	51 32 321	51 32 501	51 32 751	51 32 001	51 32 131
blau / blue	10	51 32 322	51 32 502	51 32 752	51 32 002	51 32 132
rot / red	10	51 32 323	51 32 503	51 32 753	51 32 003	51 32 133
grün / green	10	51 32 324	51 32 504	51 32 754	51 32 004	51 32 134
gelb / yellow	10	51 32 325	51 32 505	51 32 755	51 32 005	51 32 135

This document was sent online and is valid without signature

Erstellt durch / Form edited by:	Aleksandra Mil / Auszubildende	Freigegeben von / Form approved by:	H. Bauer / QS/Technik
-------------------------------------	--------------------------------	--	-----------------------