
	Qualitätsmanagement Handbuch	Formular P Form P	Seite page Seite 1 von 1
	Produkt-Spezifikation Product Specification	Änd. Datum Date of change 2022-12-05	Änd. Stand Revision 4

Allgemeine Produktdaten

96-Well PCR-Racks, Arbeits- und Aufbewahrungsstation für 0,2 und 0,5 ml PCR-Einzelgefäße, PCR-Streifen und PCR-Segmente, 8 x 12 Stellplätze im Mikrottestplattenformat, autoklavierbar, lieferbar in fünf Farben zur einfachen Kategorisierung von verschiedenen Probetypen
Transparenter Deckel in zwei Höhen lieferbar.



General Product Information

96-Well PCR-Racks, work and storage rack for 0.2 and 0.5 ml PCR tubes, PCR strips, and PCR segments, 8 x 12 positions, micro test plate format, autoclavable, available in five colors for the convenient identification of different types of samples
Transparent lids available in two different heights.

Material / material: Polypropylen / Polypropylene (PP)

Material Deckel / Material of cap: Polyethylen, Styrol-Acrylnitril / Polyethylene, Styrene-Acrylonitrile

Warentarif-Nummer / Intrastat Number- HS Code: 39269097

Sterilisiert / pre-sterilized: Ja/Yes Nein/No

CE-Zeichnung / CE-marking: Ja/Yes Nein/No

Artikel Product	Material Material	Verpackungseinheit Quantity per Pack	Bestell-Nr. Order No.
PCR-Rack weiß / white	PP	10	61 14 100
PCR-Rack gelb / yellow	PP	10	61 14 101
PCR-Rack blau / blue	PP	10	61 14 102
PCR-Rack grün / green	PP	10	61 14 103
PCR-Rack rot / red	PP	10	61 14 104
Mix-Pack aus 2x5 Farben	PP	10	61 14 105
Deckel für / lid for PCR-Rack, Höhe / height 14 mm	PP	10	61 14 110
Deckel für / lid for PCR-Rack, Höhe / height 28 mm	SAN	10	61 14 111

This document was sent online and is valid without signature

Erstellt durch / Form edited by:	Aleksandra Mil / Auszubildende	Freigegeben von / Form approved by:	H. Bauer / QS/Technik
-------------------------------------	--------------------------------	--	-----------------------