

 <small>disposables for sciences</small>	Qualitätsmanagement Handbuch	Formular P Form P	Seite page Seite 1 von 1
	Produkt-Spezifikation Product Specification	Änd. Datum Date of change 2022-12-05	Änd. Stand Revision 4

Allgemeine Produktdaten

Serologische Pipetten: RNase-, DNase-, Pyrogen-frei, TSE-, BSE-frei, einzeln steril verpackt, permanent aufgedruckte schwarze Graduierung, mit farbcodiertem Ring zur einfachen Identifikation, Farben nach Norm-Code, zusätzliche Kapazität durch negative Graduierung. Genauigkeit +/- 2-3% auf volle Volumen.

General Product Information

Serological Pipets: free from pyrogen, RNase, DNase, TSE/BSE, sterile, individually packed, permanent black-embossed graduations, with color-coded ring for convenient identification, standard color code, negative graduation for additional capacity. Accuracy is +/- 2-3%.

Material / material: Polystyrol / polystyrene (PS)

Warentarif-Nummer / Intrastat Number- HS Code: 39269097

CE-Kennzeichnung / CE-marking: Ja/Yes Nein/No

Sterilisiert / pre-sterilized: Ja/Yes Nein/No

Positives Volumen ml / Positive Volume ml	Negatives Volumen mL / Negative Volume mL	Graduierung mL / Graduation mL	Länge mm/ Length mm	Verpackungsfo rm / Packaging	Bestellnr. / Order No.
1,0 ml	0.04	0.01	279,40	1000	31 010 11
2,0 ml	0.08	0.01	282,58	1000	31 020 11
5,0 ml	2.5	0.1	348,15	200	31 050 11
10 ml	3.0	0.1	346,81	200	31 100 11
25 ml	7.0	0.25	310,90	200	31 250 11
50 ml	10.0	1.0	368,17	100	31 500 11

This document was sent online and is valid without signature

Erstellt durch / Form edited by:	Aleksandra Mil / Auszubildende	Freigegeben von / Form approved by:	H. Bauer / QS/Technik
-------------------------------------	--------------------------------	--	-----------------------