

 <small>disposables for sciences</small>	Qualitätsmanagement Handbuch	Formular P Form P	Seite page Seite 1 von 1
	Produkt-Spezifikation Product Specification	Änd. Datum Date of change 2022-12-05	Änd. Stand Revision 4

## Allgemeine Produktdaten

**ratiolab® thermocult Brutschrank**, kompakter Brutschrank für die Bakteriologie mit fest eingestellter Temperatur  $35 \pm 2^\circ\text{C}$  und Kontrollthermometer, Glastür für visuelle Kontrolle, auf Tischen aufstellbar, an Wänden und unter Schränken aufhängbar, für bis zu 20 Objektkulturen oder 8 Petrischalen,  
 Nutzraum (B x T x H): 210 x 145 x 110 mm,  
 Maße außen (B x T x H): 300 x 170 x 150 mm,  
 Gewicht: ca. 1,3 Kg,  
 Netzanschluss: 220V, 30W

## General Product Information

**ratiolab® thermocult Incubator**, Compact incubator for bacteriology with fixed temperature setting at  $35 \pm 2^\circ\text{C}$  and built-in control thermometer, glass-clear door for visual control, the incubator can be placed on work benches or be mounted on walls or underneath cabinets, for up to 20 tissue cultures or 8 petri dishes,  
 inside size (W x D x H): 210 x 145 x 110 mm,  
 outside size (W x D x H): 300 x 170 x 150 mm  
 weight: approx. 1.3 Kg,  
 power: 220V, 30W

Warentarif-Nummer / Intrastat Number- HS Code: 85143080

Sterilisiert / pre-sterilized: Ja/Yes  Nein/No

CE-Zeichnung / CE-marking: Ja/Yes  Nein/No

Artikel Product	Verpackungseinheit Quantity per Pack	Bestell-Nr. Order No.
ratiolab® thermocult	1	90 22 000
Einlegeboden / shelf	1	90 22 001

This document was sent online and is valid without signature

Erstellt durch / Form edited by:	Aleksandra Mil / Auszubildende	Freigegeben von / Form approved by:	H. Bauer / QS/Technik
-------------------------------------	--------------------------------	--	-----------------------