

 <small>disposables for sciences</small>	Qualitätsmanagement Handbuch	Formular P <i>Form P</i>	Seite <i>page</i> Seite 1 von 1
	Produkt-Spezifikation <i>Product Specification</i>	Änd. Datum <i>Date of change</i> 2025.01.21	Änd. Stand <i>Revision</i> 5

Allgemeine Produktdaten

96-Well- und 384-Well PCR-Platten

Optimaler Wärmetransfer bei kurzen Zykluszeiten, die glatten Innenseiten der Kavitäten minimieren die Bindung von Enzymen und Nukleinsäuren, zur Verhinderung von Cross-Kontamination, mit oder ohne Rand, alphanumerische Codierung zur optimalen Probenidentifikation. DNA, DNase- und RNase-frei.

General Product Information

96-well- and 384-well PCR Plates

optimal heat transfer in short cycle times, the smooth inside walls of the cavities minimize the binding of enzymes and nucleic acids, without or with edges, to minimize cross-contamination, alphanumeric coding for optimal sample identification. The plates are free of DNA, DNase and RNase.

Material / *Material*: Polypropylen / *Polypropylene*

Warentarif-Nummer / *Intrastat Number- HS Code*: 39269097

Sterilisiert / *pre-sterilized*: Ja/Yes Nein/No

CE-Zeichnung / *CE-marking*: Ja/Yes Nein/No

Volumen <i>Volume</i>	Ausführung <i>Description</i>	Packungseinheit <i>Packaging</i>	Verpackungseinheit <i>Quantity per Unit</i>	Bestell Nr. <i>Order No.</i>
96 x 0,2 ml	Mit Rand, hohes Profil / <i>with rim, high form</i>	Beutel / <i>bags</i> 5 x 10	50	86 10 011
96 x 0,2 ml	Halbrand/Semi-skirted, straight side	Beutel / <i>bags</i> 5 x 10	50	86 10 012
96 x 0,1 ml	Halbrand, FAST Type, niedriges Profil / semi-skirted, FAST Type, low profile	Beutel / <i>bags</i> 5 x 10	50	86 10 013
96 x 0,2 ml	Ohne Rand, niedriges Profil / <i>without rim, low form</i>	Beutel / <i>bags</i> 5 x 10	50	86 10 014
96 x 0,2 ml	Halbrand erhöht, ABI-Type / semi skirted, ABI Type, raised rim	Beutel / <i>bags</i> 5 x 10	50	86 10 016
96 x 0,2 ml	Halbrand, Light Cyclor Type, niedriges Profil / semi-skirted, Light Cyclor Type, low profile	Beutel / <i>bags</i> 5 x 10	50	86 10 017
384 x 0,04 ml	Vollrand / Fully skirted	Beutel / <i>bags</i> 5 x 10	50	86 10 020
384 x 0,04 ml	Vollrand, weiß, schwarz kodiert / Fully skirted, white, black coded	Beutel / <i>bags</i> 5 x 10	50	86 10 021

This document was sent online and is valid without signature.

Erstellt durch / <i>Form edited by:</i>	E. Tagik QS/ Technik	Freigegeben von / <i>Form approved by:</i>	E. Tolmach/ Einkauf
--------------------------------------------	----------------------	-----------------------------------------------	---------------------